



DÜRR DENTAL SE · Postfach 1264 · 74302 Bietigheim-Bissingen

TO WHOM IT MAY CONCERN

DÜRR DENTAL SE
Höpfigheimer Str. 17
74321 Bietigheim-Bissingen
Deutschland / Germany
Tel +49 7142 705-0
Fax +49 7142 705-500
E-Mail info@duerrdental.com
Internet www.duerrdental.com

Unsere Zeichen / Our Ref. Tel Fax E-Mail

Datum / Date
21 April 2020

Statement on the virucidal activity of products from Dürr System-Hygiene against coronavirus (SARS-CoV-2) according to the European standards EN 14476 and EN 17111

Coronaviruses belong to the group of enveloped viruses which are generally easy to inactivate. In Europe, the disinfectant has to be tested against the specified test viruses (vacciniavirus or adenovirus/murine norovirus or adenovirus/murine norovirus/poliovirus, respectively) according to EN 14476 and EN 17111 (additionally required in case of instrument disinfection) to claim the virucidal activity against all enveloped viruses. According to EN 14476 and EN 17111 these test viruses serve as surrogate viruses since several viruses, such as HBV, HCV or HIV are not culturable in common cell cultures. Therefore, the usual disinfectant tests cannot be performed with these viruses.

Therefore, after successful experiments with the above mentioned non-enveloped viruses or enveloped virus, the disinfectant is also effective against all enveloped viruses including HBV, HCV and HIV (blood-borne viruses) as well as against coronaviruses.

We hereby confirm that the following products from Dürr System-Hygiene fulfill the European virucidal standards (EN 14476 and EN 17111, respectively) and are effective against all enveloped viruses as well as against the new type of coronavirus (SARS-CoV-2).

Instruments:

- ID 212/ID 212 forte Instrument disinfection
- ID 213 Instrument disinfection
- ID 220 Bur disinfection

Surfaces:

- FD 300/FD 300 top wipes Surface disinfection
- FD 312 Surface disinfection
- FD 322/FD 322 premium wipes/FD 322 top wipes Quick-acting disinfection
- FD 333/FD 333 wipes/FD 333 top wipes Quick-acting disinfection
- FD 333 forte/FD 333 forte premium wipes/FD 333 forte wipes Quick-acting disinfection
- FD 350 Disinfection wipes
- FD 366 sensitive/FD 366 sensitive wipes/FD 366 sensitive top wipes Disinfection of sensitive surfaces

Skin and hands:

- HD 410 Hand disinfection
- HD 412 essential Hand disinfection

Special areas:

- Orotol® plus Disinfection of suction systems
- Orotol® ultra Disinfection of suction systems
- MD 520 Impression disinfection

A handwritten signature in blue ink, appearing to read 'Pflug'.

Christian Pflug
Vice President
Hygiene and Dental Care

A handwritten signature in blue ink, appearing to read 'i.v. M. Diederich'.

Marc Diederich
Head of Product Management
Hygiene and Dental Care